

Product datasheet for **MP216374**

Sult2b1 Mouse qPCR Primer Pair (NM_017465)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54200
Forward Sequence:	GCAGTTTGCTCCTGGTTTGAC
Reverse Sequence:	CCAGGAATCACAGATGCGTTG
ACCN:	<u>BC009811</u> , <u>BC009813</u> , <u>NM_017465</u> , <u>NM_017465.1</u> , <u>NM_017465.2</u> , <u>BB867946</u> , <u>BM899747</u> , <u>NM_017465.3</u>
UniProt ID:	<u>O35400</u>
Synonyms:	AI326997; BB173635; ST2B1; SULT2B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)