

Product datasheet for **MP216385**

Sod1 Mouse qPCR Primer Pair (NM_011434)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20655
Forward Sequence:	GGTGAACCAAGTTGTGTTGTCAGG
Reverse Sequence:	ATGAGGTCCTGCACTGGTACAG
ACCN:	BC002066 , BC048874 , BC086886 , NM_011434 , NM_011434.1 , BC032986 , BC057074 , BC069974 , NM_011434.2
UniProt ID:	P08228
Synonyms:	B430204E11Rik; Cu/Zn-SOD; CuZnSOD; lpo-1; lpo1; Sod-1; SODC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)