

Product datasheet for **MP216415**

Suv39h2 Mouse qPCR Primer Pair (NM_022724)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64707
Forward Sequence:	CCATAAACGCTGGAGAAGAGCTG
Reverse Sequence:	CTGCAAGTCTCGGCTCCACATT
ACCN:	<u>BC032960, NM_022724, NM_022724.1, NM_022724.2, NM_022724.3, NM_022724.4, BY752501</u>
UniProt ID:	<u>Q9EQQ0</u>
Synonyms:	4930507K23Rik; AA536750; D030054H19Rik; D2Ertd544e; KMT1B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)