

Product datasheet for **MP216444**

Smarca4 Mouse qPCR Primer Pair (NM_011417)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20586
Forward Sequence:	GAAAGTGGCTCTGAAGAGGAGG
Reverse Sequence:	TCCACCTCAGAGACATCATCGC
ACCN:	BC079560 , NM_011417 , NM_011417.1 , NM_011417.2 , NM_011417.3 , BC014760 , BC023186 , BC023830 , BC026672 , BC032162 , BC054056 , BC060229 , BC061214
UniProt ID:	Q3TKT4
Synonyms:	b2b508.1Clo; b2b692Clo; Brg1; HP1-BP72; SNF2beta; SW1/SNF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)