

Product datasheet for **MP216445**

Smarca5 Mouse qPCR Primer Pair (NM_053124)

Product data:

| | |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 93762 |
| Forward Sequence: | GCTGCTCAGAAGACTCCAACCT |
| Reverse Sequence: | TCCTCTTGCTCCGTTCTACGGT |
| ACCN: | <u>BC053069</u> , <u>NM_053124</u> , <u>NM_053124.1</u> , <u>NM_053124.2</u> , <u>BC009130</u> , <u>BC021922</u> |
| UniProt ID: | <u>Q91ZW3</u> |
| Synonyms: | 4933427E24Rik; D030040M08Rik; D330027N15Rik; MommeD4; Snf2h |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)