

Product datasheet for **MP216513**

Syvn1 Mouse qPCR Primer Pair (NM_028769)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74126
Forward Sequence:	CCAACATCTCCTGGCTCTTCCA
Reverse Sequence:	CAGGATGCTGTGATAAGCGTGG
ACCN:	<u>BC042199</u> , <u>BC046829</u> , <u>BC057917</u> , <u>BC080722</u> , <u>NM_028769</u> , <u>NM_028769.1</u> , <u>NM_028769.2</u> , <u>NM_028769.3</u> , <u>NM_028769.4</u> , <u>NM_028769.5</u> , <u>BB676375</u> , <u>BC032931</u>
UniProt ID:	<u>Q9DBY1</u>
Synonyms:	1200010C09Rik; AW211966; C85322; D530017H19Rik; Hrd1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)