

Product datasheet for **MP216671**

Tbc1d24 Mouse qPCR Primer Pair (NM_173186)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	224617
Forward Sequence:	CTCCTGCTCTTCTCATCACTGC
Reverse Sequence:	CCAGTCTGTTGACAGGTAAGCC
ACCN:	<u>BC080845, BC094417, NM_173186, NM_173186.1, NM_173186.2, NM_173186.3, NM_173186.4</u>
UniProt ID:	<u>Q3UUG6</u>
Synonyms:	9630033P11; C530046L02Rik; mKIAA1171
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)