

Product datasheet for **MP216702**

Timd2 Mouse qPCR Primer Pair (NM_134249)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	171284
Forward Sequence:	CTCTACTTCTCCAACACCAGCAC
Reverse Sequence:	CCAAGGGTCATCTGAGGATGTC
ACCN:	<u>BC025096</u> , <u>BC028829</u> , <u>NM_134249</u> , <u>NM_134249.1</u> , <u>NM_134249.2</u> , <u>NM_134249.3</u> , <u>NM_134249.4</u> , <u>NM_134249.5</u> , <u>BY731803</u>
UniProt ID:	<u>Q8R183</u>
Synonyms:	AA816106; AU019358; AU019662; C78111; G9D3; TIM-2; Tim2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)