

Product datasheet for **MP216723**

Tcp10a Mouse qPCR Primer Pair (NM_011553)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21461
Forward Sequence:	GAGGAACAGACTCCCAGGATGA
Reverse Sequence:	TCAGCCAGTGACTCTGTGTCCA
ACCN:	<u>BC050771</u> , <u>BC052353</u> , <u>NM_011553</u> , <u>NM_011553.1</u> , <u>NM_011553.2</u> , <u>NM_011553.3</u> , <u>NM_011553.4</u> , <u>BC152833</u>
Synonyms:	D17Leh66A; D17Leh66aa; Gm10326; T66A-a; Tcp-10a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)