

## Product datasheet for **MP216738**

### Tead1 Mouse qPCR Primer Pair (NM\_009346)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21676
Forward Sequence:	TCCGCTTCCTTGAACAGCAGAG
Reverse Sequence:	GGTCACTGTAAGAAATGGTTGGC
ACCN:	<u><a href="#">BC060138</a></u> , <u><a href="#">NM_009346</a></u> , <u><a href="#">NM_009346.1</a></u> , <u><a href="#">NM_009346.3</a></u> , <u><a href="#">BC027158</a></u> , <u><a href="#">BC049153</a></u> , <u><a href="#">BC052021</a></u> , <u><a href="#">BG064848</a></u> , <u><a href="#">BM214769</a></u>
Synonyms:	2610024B07Rik; B230114H05Rik; Gtrgeo5; mTEF-1; Tcf13; TEAD-1; TEF-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)