

Product datasheet for **MP216821**

Tspan14 Mouse qPCR Primer Pair (NM_145928)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 52588 |
| Forward Sequence: | GCAGAACCTCATTGACTCCCTTC |
| Reverse Sequence: | ATCTGGCACACAGCAGGAGAAG |
| ACCN: | BC024611 , BC025568 , BC026574 , NM_145928 , NM_145928.1 , NM_145928.2 , BG807867 , BY209307 |
| UniProt ID: | Q8QZY6 |
| Synonyms: | AA958793; D14Ertd226e; Tm4sf14; tspan-14 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)