

Product datasheet for **MP216831**

Tspan5 Mouse qPCR Primer Pair (NM_019571)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56224
Forward Sequence:	GCACAGATTCCAATGCAAGCCG
Reverse Sequence:	GGTTTCTGCCTGGCATCATAGC
ACCN:	<u>BC046784, BC058695, NM_019571, NM_019571.1, NM_019571.2, NM_019571.3, NM_019571.4, NM_019571.5</u>
UniProt ID:	<u>P62080</u>
Synonyms:	2810455A09Rik; 4930505M03Rik; NET-4; Tm4sf9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)