

Product datasheet for **MP216874**

Thap11 Mouse qPCR Primer Pair (NM_021513)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 59016 |
| Forward Sequence: | GCTGTGCTTCTTACGCTTCAGG |
| Reverse Sequence: | ACTCGACTTGCACTGTGAGGTC |
| ACCN: | <u>NM_021513</u> , <u>NM_021513.1</u> , <u>NM_021513.2</u> , <u>BC138966</u> , <u>BC003949</u> , <u>BC138987</u> |
| UniProt ID: | <u>Q9JJD0</u> |
| Synonyms: | 2810036E22Rik; AB041579; CTG-B43a; CTG-B45d; Ronin |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)