

Product datasheet for **MP217011**

Tnip1 Mouse qPCR Primer Pair (NM_021327)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57783
Forward Sequence:	GCTGGAAGTGAACAAGCAGTGG
Reverse Sequence:	CCAGCTCAGTTACCTGTTTCTGC
ACCN:	<u>BC008186</u> , <u>BC008665</u> , <u>NM_021327</u> , <u>NM_021327.1</u> , <u>NM_021327.2</u> , <u>NM_021327.4</u> , <u>BQ030862</u> , <u>BY170310</u> , <u>BY197901</u>
UniProt ID:	<u>Q9WUU8</u>
Synonyms:	ABIN; ABIN-1; ABIN1; AU018810; Naf1; Nef; VAN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)