

Product datasheet for **MP217018**

Tirap Mouse qPCR Primer Pair (NM_054096)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	117149
Forward Sequence:	ATCTCCCAGGAAAGCCACCTCT
Reverse Sequence:	GGTAGGTGACATTCTGAACTGC
ACCN:	BC005584 , NM_054096 , NM_054096.1 , NM_054096.2 , BY507516 , BY740695 , BY743490
UniProt ID:	Q99JY1
Synonyms:	AA407980; C130027E04Rik; Mal; Tlr4ap; Wyatt
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)