

Product datasheet for **MP217031**

Ticam1 Mouse qPCR Primer Pair (NM_174989)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	106759
Forward Sequence:	ATCCATGCCAGGGCTGATGAAC
Reverse Sequence:	CGATGGCATCTTGGAGACAGTG
ACCN:	BC094338 , NM_174989 , NM_174989.1 , NM_174989.2 , NM_174989.3 , NM_174989.4 , BC033406 , BC037048 , BC062191 , NM_174989.5
UniProt ID:	Q80UF7
Synonyms:	AW046014; AW547018; TICAM-1; TRIF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)