

Product datasheet for **MP217038**

Topbp1 Mouse qPCR Primer Pair (NM_176979)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	235559
Forward Sequence:	CGGAACTAAGCTCACCCTCTC
Reverse Sequence:	CAACGCTGTCAAGTGGTTCCG
ACCN:	<u>BC049797</u> , <u>NM_176979</u> , <u>NM_176979.3</u> , <u>NM_176979.4</u> , <u>NM_176979.5</u> , <u>BC006707</u> , <u>BC007170</u> , <u>BC026608</u> , <u>BC062111</u>
UniProt ID:	<u>Q6ZQF0</u>
Synonyms:	1110031N14Rik; 2810429C13Rik; AI256758; D430026L04Rik; mKIAA0259
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)