

## Product datasheet for **MP217048**

### Tlk1 Mouse qPCR Primer Pair (NM\_172664)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	228012
Forward Sequence:	GGAACCTTCACATCCGAGAGCT
Reverse Sequence:	TAGACCTCACTAAAGCCGCCTC
ACCN:	<a href="#">BC051641</a> , <a href="#">NM_172664</a> , <a href="#">NM_172664.1</a> , <a href="#">NM_172664.2</a> , <a href="#">NM_172664.3</a> , <a href="#">BC033931</a> , <a href="#">BC047088</a> , <a href="#">BC057369</a>
UniProt ID:	<a href="#">Q8C0V0</a>
Synonyms:	4930545J15Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)