

Product datasheet for **MP217120**

Tcf4 Mouse qPCR Primer Pair (NM_013685)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21413
Forward Sequence:	CCTCCAATCCTTCAACTCCTGTG
Reverse Sequence:	TCCAAACGGTCTTCGATTCGGC
ACCN:	<u>BC014293</u> , <u>BC043050</u> , <u>NM_013685</u> , <u>NM_013685.1</u> , <u>NM_013685.2</u> , <u>BC094350</u> , <u>BY770794</u>
UniProt ID:	<u>Q60722</u>
Synonyms:	5730422P05Rik; ASP-I2; bHLHb19; E2-2; E2.2; ITF-2; ITF-2b; ITF2; ME2; MITF-2A; MITF-2B; SEF-2; SEF2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)