

Product datasheet for **MP217126**

Tfap2d Mouse qPCR Primer Pair (NM_153154)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226896
Forward Sequence:	CTATCGCTGAGGTAAGAGGCG
Reverse Sequence:	CTCTCAGGCAACGACCACCATT
ACCN:	<u>NM_153154</u> , <u>NM_153154.1</u> , <u>NM_153154.2</u> , <u>NM_153154.3</u> , <u>BC152324</u> , <u>BB676311</u> , <u>BC152325</u> , <u>BC160252</u>
UniProt ID:	<u>Q91ZK0</u>
Synonyms:	Tcfap2d
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)