

Product datasheet for **MP217136**

Tcf3 Mouse qPCR Primer Pair (NM_011548)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21423
Forward Sequence:	CCATGCTAGGTGACGGCTCTTC
Reverse Sequence:	GCGAGCCATTAACCTCAGATCC
ACCN:	<u>BC006860, BC018260, NM_011548, NM_011548.1, NM_011548.2, NM_011548.3, NM_011548.4</u>
UniProt ID:	<u>P15806</u>
Synonyms:	A1; AA408400; ALF2; AW209082; bHLHb21; E2A; E12; E12/E47; E47; KA1; ME2; Pan1; Pan2; TCF-3; Tcf2a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)