

Product datasheet for **MP217137**

Tfe3 Mouse qPCR Primer Pair (NM_172472)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	209446
Forward Sequence:	CTATCTTCAGGAGGCACTGCA
Reverse Sequence:	CTCGCGTTTGATGTTAGGCAGC
ACCN:	BC056358 , BC063047 , NM_172472 , NM_172472.1 , NM_172472.2 , NM_172472.3 , BC032963 , BC051568 , BY208269 , BY299767
UniProt ID:	Q64092
Synonyms:	bHLHe33; F830016E06Rik; Tcf3; Tfe-3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)