

## Product datasheet for **MP217148**

### Sin3b Mouse qPCR Primer Pair (NM\_009188)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Gene ID:          | 20467                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | TACGAGGAGCAACTGCATCGCT                                                                                                                                                                                                                                                                                                                        |
| Reverse Sequence: | TCAGGAGCCATCCGTGACAACT                                                                                                                                                                                                                                                                                                                        |
| ACCN:             | <a href="#">BC129859</a> , <a href="#">NM_009188</a> , <a href="#">NM_009188.1</a> , <a href="#">NM_009188.2</a> , <a href="#">NM_009188.3</a> , <a href="#">NM_009188.4</a> , <a href="#">BC020049</a> , <a href="#">BC021160</a> , <a href="#">BC051536</a> , <a href="#">BC096367</a> , <a href="#">BF318403</a>                           |
| UniProt ID:       | <a href="#">Q62141</a>                                                                                                                                                                                                                                                                                                                        |
| Synonyms:         | 2810430C10Rik                                                                                                                                                                                                                                                                                                                                 |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |

[View online »](#)