

## Product datasheet for **MP217178**

### Mdm4 Mouse qPCR Primer Pair (NM\_008575)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17248
Forward Sequence:	GCTCTCGCACAGGATCACACTA
Reverse Sequence:	ATGTCGTGAGGTAGGCAGTGTG
ACCN:	<a href="#">BC003750</a> , <a href="#">NM_008575</a> , <a href="#">NM_008575.1</a> , <a href="#">NM_008575.2</a> , <a href="#">NM_008575.3</a> , <a href="#">NM_008575.4</a> , <a href="#">BC023912</a> , <a href="#">BC025868</a> , <a href="#">BC058680</a>
UniProt ID:	<a href="#">O35618</a>
Synonyms:	4933417N07Rik; AA414968; AL023055; AU018793; AU021806; C85810; Mdmx
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)