

## Product datasheet for **MP217179**

### Tra2a Mouse qPCR Primer Pair (NM\_198102)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 101214   |
| Forward Sequence: | AACGGGAAGCTCTGCTCGTGTA   |
| Reverse Sequence: | GTGTAGCGTCTATGCGAATGCC   |
| ACCN:             | <u><a href="#">BC057448</a></u> , <u><a href="#">NM_198102</a></u> , <u><a href="#">NM_198102.1</a></u> , <u><a href="#">NM_198102.2</a></u> , <u><a href="#">BB873436</a></u> , <u><a href="#">BC023018</a></u> , <u><a href="#">BC058764</a></u> , <u><a href="#">BE134105</a></u>   |
| Synonyms:         | 1500010G04Rik; AL022798; G430041M01Rik; mAWMS1   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)