

Product datasheet for **MP217209**

Trpc3 Mouse qPCR Primer Pair (NM_019510)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22065
Forward Sequence:	ATTCCTCGCCATCGGCTATTGG
Reverse Sequence:	GAGGCGTTGAATACAAGCAGACC
ACCN:	<u>BC108989</u> , <u>NM_019510</u> , <u>NM_019510.1</u> , <u>NM_019510.2</u> , <u>BC108988</u>
Synonyms:	Mwk; Trcp3; Trp3; Trpp3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)