

Product datasheet for **MP217217**

Trpm3 Mouse qPCR Primer Pair (NM_001035239)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226025
Forward Sequence:	TCTACAACGCACCCATCGTGAAG
Reverse Sequence:	GCGCTCCATCTTACTAACACG
ACCN:	<u>NM_001035239</u> , <u>NM_001035239.1</u> , <u>NM_001035239.2</u> , <u>BC141400</u> , <u>BB467120</u> , <u>BC037449</u> , <u>BC141412</u> , <u>BE995087</u> , <u>NM_001035239.3</u>
Synonyms:	6330504P12Rik; 9330180E14; AU018608; B930001P07Rik; LTRPC3; MLSN2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)