

Product datasheet for **MP217279**

Tmco3 Mouse qPCR Primer Pair (NM_172282)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	234076
Forward Sequence:	CACATGGTCACCGAGGAGATCA
Reverse Sequence:	TGAGGAACACCAGCACCGTCAA
ACCN:	<u>BC039955</u> , <u>NM_172282</u> , <u>NM_172282.1</u> , <u>NM_172282.2</u> , <u>BC026890</u> , <u>BC053487</u> , <u>BY146007</u>
UniProt ID:	<u>Q8BH01</u>
Synonyms:	B230339H12Rik; C87304
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)