

Product datasheet for **MP217398**

Ano6 Mouse qPCR Primer Pair (NM_175344)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	105722
Forward Sequence:	CGCATCCCATTTACCACCTGTG
Reverse Sequence:	CGATGAACACGGACAGCCTGTA
ACCN:	BC060732 , NM_175344 , NM_175344.1 , NM_175344.2 , NM_175344.3 , NM_175344.4 , BC023745 , BC052511 , BE853067 , BY717873
UniProt ID:	Q6P9J9
Synonyms:	2900059G15Rik; AA407480; AW554778; F730003B03Rik; Tmem16f
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)