

Product datasheet for **MP217525**

Tcof1 Mouse qPCR Primer Pair (NM_011552)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21453
Forward Sequence:	GAGGACTCAGACAGCAGTAGTG
Reverse Sequence:	GTGTTACAGCCTGGACTCCTT
ACCN:	BC021514 , BC057342 , NM_011552 , NM_011552.1 , NM_011552.2 , NM_011552.3 , BC021514.1 , BC046608 , BC052669 , BC060105
UniProt ID:	O08784
Synonyms:	AA408847; AW209012; treacle
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)