

## Product datasheet for **MP217538**

### Trem2 Mouse qPCR Primer Pair (NM\_031254)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83433
Forward Sequence:	CTACCAAGTGTCTCAGAGTCTCCGA
Reverse Sequence:	CCTCGAAACTCGATGACTCCTC
ACCN:	<a href="#">BC032959</a> , <a href="#">BC033485</a> , <a href="#">BC052784</a> , <a href="#">NM_031254</a> , <a href="#">NM_031254.1</a> , <a href="#">NM_031254.2</a> , <a href="#">NM_031254.3</a> , <a href="#">BY541740</a>
UniProt ID:	<a href="#">Q99NH8</a>
Synonyms:	Trem; TREM-2; Trem2a; Trem2b; Trem2c
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)