

Product datasheet for **MP217563**

Trim2 Mouse qPCR Primer Pair (NM_030706)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80890
Forward Sequence:	CACGAACCTCATGGATGTGCTG
Reverse Sequence:	GGTCTCACAGGACTGGCAGTAA
ACCN:	<u>BC058961, NM_030706, NM_030706.1, NM_030706.2, NM_030706.3, BC058668, BY124765</u>
UniProt ID:	<u>Q9ESN6</u>
Synonyms:	mKIAA0517; narf
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)