

Product datasheet for **MP217567**

Trim25 Mouse qPCR Primer Pair (NM_009546)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	217069
Forward Sequence:	ACAAGTTGCCACCTTTGGAGC
Reverse Sequence:	CACCTTAGCCTTGACTCCTACAG
ACCN:	BC006908 , BC034276 , NM_009546 , NM_009546.1 , NM_009546.2 , BC028430
UniProt ID:	Q61510
Synonyms:	AA960166; AL022677; EFP; Zfp147
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)