

Product datasheet for **MP217638**

Tpm1 Mouse qPCR Primer Pair (NM_024427)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22003
Forward Sequence:	TGAGCAAGCGGAGGCTGATAAG
Reverse Sequence:	GCATCTTTGAGAGCCTCGGAGT
ACCN:	BC132035 , BC132037 , NM_024427 , NM_024427.1 , NM_024427.2 , NM_024427.3 , NM_024427.4 , BC145307 , BC026720 , BC131643
UniProt ID:	P58771
Synonyms:	AA986836; AI854628; alpha-TM; TM2; Tm3; Tmpa; Tpm-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)