

Product datasheet for **MP217659**

Tph1 Mouse qPCR Primer Pair (NM_009414)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21990
Forward Sequence:	TGTTGACTGCGACATCAGCCGA
Reverse Sequence:	GGAAACCAAGGGACAGTCTCCA
ACCN:	<u>BC072582</u> , <u>NM_009414</u> , <u>NM_009414.1</u> , <u>NM_009414.2</u> , <u>NM_009414.3</u> , <u>BC038864</u> , <u>BM932391</u>
UniProt ID:	<u>P17532</u>
Synonyms:	Tph
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)