

Product datasheet for **MP217667**

Tsc22d4 Mouse qPCR Primer Pair (NM_023910)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	78829
Forward Sequence:	AGAGATGGAGCAGTTCGGCTGA
Reverse Sequence:	TATGGACCAGGCTGGCATCGAA
ACCN:	<u>BC018544</u> , <u>BC023761</u> , <u>NM_023910</u> , <u>NM_023910.1</u> , <u>NM_023910.2</u> , <u>NM_023910.3</u> , <u>NM_023910.4</u> , <u>NM_023910.5</u> , <u>NM_023910.6</u> , <u>BC145721</u> , <u>BC145723</u> , <u>BE685784</u> , <u>BY750078</u>
UniProt ID:	<u>Q9EQN3</u>
Synonyms:	0610009M14Rik; 1700023B23Rik; AI415410; Thg-1pit; Tilz2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)