

Product datasheet for **MP217668**

Tsc22d1 Mouse qPCR Primer Pair (NM_207652)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21807
Forward Sequence:	AAGGTGCTACCGCTGACAACAC
Reverse Sequence:	CCTCACCGCATACATCAAATGGC
ACCN:	<u>BC058660</u> , <u>NM_207652</u> , <u>NM_207652.1</u> , <u>NM_207652.2</u> , <u>NM_207652.3</u> , <u>BC026703</u> , <u>BY440762</u> , <u>BY720238</u>
UniProt ID:	<u>P62500</u>
Synonyms:	AA589566; AW105905; Egr5; Tgfb1i4; Tsc; TSC-22; Tsc22
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)