

Product datasheet for **MP217674**

Tsku Mouse qPCR Primer Pair (NM_001024619)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	244152
Forward Sequence:	CACTTGACCTCAGCCACAATGG
Reverse Sequence:	GACCTCTCGTAGTCGGTTATGG
ACCN:	<u>BC116310</u> , <u>BC116311</u> , <u>BC129977</u> , <u>BC129978</u> , <u>NM_001024619</u> , <u>NM_001024619.1</u> , <u>NM_001024619.2</u> , <u>NM_001024619.3</u>
UniProt ID:	<u>Q8CBR6</u>
Synonyms:	9530051K01Rik; E2ig4; Lrrc54; Tsk
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)