

Product datasheet for **MP217676**

Tufm Mouse qPCR Primer Pair (NM_172745)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	233870
Forward Sequence:	TGGATGGCTGTATCCTGGTGGT
Reverse Sequence:	CACCATCTCTGAGTCCTGGACA
ACCN:	<u>BC060959</u> , <u>BC100596</u> , <u>NM_172745</u> , <u>NM_172745.1</u> , <u>NM_172745.2</u> , <u>NM_172745.3</u> , <u>BU699735</u> , <u>BY246546</u>
UniProt ID:	<u>Q8BFR5</u>
Synonyms:	2300002G02Rik; C76308; C76389; EF-TuMT; EFTU
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)