

Product datasheet for **MP217731**

Tdrd3 Mouse qPCR Primer Pair (NM_172605)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	219249
Forward Sequence:	CTACGGCAACTATGAAGAGGTGC
Reverse Sequence:	GTAGACCTTCGTGGCTGTCCAT
ACCN:	BC005670 , BC066144 , NM_172605 , NM_172605.1 , NM_172605.2 , NM_172605.3 , BB126096 , BB668608 , BE983697
UniProt ID:	Q91W18
Synonyms:	4732418C03Rik; 6720468N18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)