

## Product datasheet for **MP217736**

### Tnfsf10 Mouse qPCR Primer Pair (NM\_009425)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 22035  |
| Forward Sequence: | GGAAGACCTCAGAAAGTGGCAG   |
| Reverse Sequence: | TTCCGAGAGGACTCCCAGGAT  |
| ACCN:             | <u><a href="#">NM_009425</a></u> , <u><a href="#">NM_009425.1</a></u> , <u><a href="#">NM_009425.2</a></u> , <u><a href="#">BC138807</a></u> , <u><a href="#">BC011469</a></u> , <u><a href="#">BC138808</a></u> , <u><a href="#">BY734631</a></u>   |
| UniProt ID:       | <u><a href="#">P50592</a></u>  |
| Synonyms:         | A330042I21Rik; AI448571; APO-2L; Ly81; TL2; Tnlg6a; Trail  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)