

Product datasheet for **MP217737**

Tnfsf11 Mouse qPCR Primer Pair (NM_011613)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21943
Forward Sequence:	GTGAAGACACACTACCTGACTCC
Reverse Sequence:	GCCACATCCAACCATGAGCCTT
ACCN:	<u>BC125603</u> , <u>BC131970</u> , <u>NM_011613</u> , <u>NM_011613.1</u> , <u>NM_011613.2</u> , <u>NM_011613.3</u>
UniProt ID:	<u>O35235</u>
Synonyms:	Ly109l; ODF; OPGL; RANKL; Trance
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)