

Product datasheet for **MP217750**

Tnfrsf10b Mouse qPCR Primer Pair (NM_020275)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21933
Forward Sequence:	TGTGTCGATGCAAACCAAGGCAC
Reverse Sequence:	GCCGTTTTGGAGACACACTTCC
ACCN:	<u>BC065141, NM_020275, NM_020275.1, NM_020275.2, NM_020275.3, NM_020275.4</u>
UniProt ID:	<u>Q9QZM4</u>
Synonyms:	DR5; KILLER; Ly98; MK; TRAILR2; TRICK2A; TRICK2B; TRICKB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)