

## Product datasheet for **MP217760**

### **Tnfrsf1a Mouse qPCR Primer Pair (NM\_011609)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	21937
<b>Forward Sequence:</b>	GTGTGGCTGTAAGGAGAACCAG
<b>Reverse Sequence:</b>	CACACGGTGTCTGAGTCTCCT
<b>ACCN:</b>	<u><a href="#">BC004599</a></u> , <u><a href="#">BC052675</a></u> , <u><a href="#">NM_011609</a></u> , <u><a href="#">NM_011609.1</a></u> , <u><a href="#">NM_011609.2</a></u> , <u><a href="#">NM_011609.3</a></u> , <u><a href="#">NM_011609.4</a></u> , <u><a href="#">BC096429</a></u>
<b>UniProt ID:</b>	<u><a href="#">P25118</a></u>
<b>Synonyms:</b>	CD120; CD120a; FPF; p5; p55; p55-R; TN; Tnf; TNF-; TNF-a; TNF-alphaR1; TNF-R; TNF-R-I; TNF-R1; TNF-R55; TNFaI; TNFalpha-R1; TNFAR; TNFR; Tnfr-2; Tnfr1; TNFR60; TNFRI; TNFRp55
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)