

Product datasheet for **MP217778**

Tpd52l1 Mouse qPCR Primer Pair (NM_009413)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21987
Forward Sequence:	GAGGAGAGGGTTGAGACAACTG
Reverse Sequence:	TGTGCTGTGGAGTTCAGGACCT
ACCN:	<u>NM_009413</u> , <u>NM_009413.1</u> , <u>BC141068</u> , <u>BC145321</u> , <u>BC145322</u> , <u>BC160227</u> , <u>BC171924</u> , <u>NM_009413.2</u>
UniProt ID:	<u>O54818</u>
Synonyms:	D53
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

