

Product datasheet for **MP217790**

Epcam Mouse qPCR Primer Pair (NM_008532)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 17075 |
| Forward Sequence: | GAGTCCGAAGAACCGACAAGGA |
| Reverse Sequence: | GATGTGAACGCCTCTTGAAGCG |
| ACCN: | <u>BC005618</u> , <u>BC094465</u> , <u>NM_008532</u> , <u>NM_008532.1</u> , <u>NM_008532.2</u> , <u>BI986346</u> |
| UniProt ID: | <u>Q99JW5</u> |
| Synonyms: | CD326; EGP; EGP-2; Egp314; Ep-CAM; EpCAM1; GA733-2; gp40; Ly74; Tacsd1; Tacstd1; TROP1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)