

Product datasheet for **MP217825**

Snrnp70 Mouse qPCR Primer Pair (NM_009224)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20637
Forward Sequence:	GCCTTCAAGACTCTGTTCGTGG
Reverse Sequence:	CTCGATGAAGGCATAACCCACGG
ACCN:	<u>BC131994</u> , <u>BC132020</u> , <u>NM_009224</u> , <u>NM_009224.2</u> , <u>NM_009224.3</u> , <u>NM_009224.4</u> , <u>NM_009224.5</u> , <u>BC002169</u> , <u>BC043938</u> , <u>BC049128</u> , <u>BC051414</u> , <u>BC058513</u> , <u>BC094067</u> , <u>BE990366</u> , <u>BF019547</u> , <u>BY297843</u>
UniProt ID:	<u>Q62376</u>
Synonyms:	2700022N21Rik; 3200002N22Rik; AI325098; R74807; Rnu; Rnu1p70; Rnulp70; Sn; Snrp70; Sr; Snrp70; U1-70; U1-70K
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)