

## Product datasheet for **MP217834**

### Methig1 Mouse qPCR Primer Pair (NM\_001024672)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                              |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                            |
| Gene ID:          | 554292                                                                                                                                                                                                                                                                                                                                       |
| Forward Sequence: | GTCTTCAGATGAGTGGTCAGCAG                                                                                                                                                                                                                                                                                                                      |
| Reverse Sequence: | GCTTGTAAGACCGTAGGACAGC                                                                                                                                                                                                                                                                                                                       |
| ACCN:             | <u><a href="#">NM_001024672</a></u> , <u><a href="#">NM_001024672.1</a></u> , <u><a href="#">NM_001024672.2</a></u> , <u><a href="#">NM_001024672.3</a></u> , <u><a href="#">BC139180</a></u> , <u><a href="#">BC139182</a></u> , <u><a href="#">BC145570</a></u>                                                                            |
| Synonyms:         | Mettl7a2-Higd1c; Mettl7a2Higd1c; Mettl7a2higd1c; UbiE-YGHL1                                                                                                                                                                                                                                                                                  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)                                                                                                                                                                                                                                               |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                              |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                 |



[View online »](#)