

Product datasheet for **MP217851**

Ubash3a Mouse qPCR Primer Pair (NM_177823)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 328795 |
| Forward Sequence: | AGAGAGTGGACCAGGTCTTTGG |
| Reverse Sequence: | CTTTGATGCCGTTGCTCCTCT |
| ACCN: | <u>NM_177823</u> , <u>NM_177823.1</u> , <u>NM_177823.2</u> , <u>NM_177823.3</u> , <u>NM_177823.4</u> , <u>BB736706</u> , <u>BC027561</u> , <u>BC140312</u> , <u>BC148683</u> , <u>BY227020</u> |
| UniProt ID: | <u>Q3V3E1</u> |
| Synonyms: | 5830413C03Rik; C330001M22; Sts-2; TULA; TULA-1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)