

Product datasheet for **MP217855**

Ubap2l Mouse qPCR Primer Pair (NM_028475)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74383
Forward Sequence:	GGTCTTCAGACAACCAGTCCTC
Reverse Sequence:	ATCTCCACAGCCAGAGCAGGAA
ACCN:	<u>BC023906</u> , <u>BC029075</u> , <u>BC050910</u> , <u>NM_028475</u> , <u>NM_028475.1</u> , <u>NM_028475.2</u> , <u>BB648932</u> , <u>BB842333</u> , <u>BM124250</u> , <u>BM935680</u> , <u>NM_028475.3</u>
UniProt ID:	<u>Q80X50</u>
Synonyms:	3110083O19Rik; 4932431F02Rik; A430103N23Rik; C77168; mKIAA0144; Nice-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)